

Digital Imaging

Alberto F. Goldszal, PhD
Acting Chief, DIN/PACS/RIS
Imaging Sciences Program, CC, NIH



DIN/PACS/RIS Definitions

DIN = Digital Imaging Network

PACS = Picture Archiving & Communication System

RIS = Radiology Information System

- PACS/RIS are medical devices, not merely information systems.
- Complex automation systems for acquiring, transmitting, managing, storing and displaying digital diagnostic images and associated text information.
- PACS/RIS systems are the application of information age technology to *improve the quality and efficiency* of radiology services.
- Enable integration of new and existing medical record data including image, biosignal, textual and demographic data toward a filmless and paperless operation.


Why PACS/RIS?

- Decrease number of repeated examinations.
- Decrease/eliminate film and chemical costs.
- Decrease turnaround time of exams.
- Improve access to prior studies.
- Improve the quality of images.
- Decrease frequency of misidentified images.
- Improve equipment availability.
- Improve efficiency of after-hours interpretation.
- Improve availability of images and reports.
- Minimize/eliminate lost films.

Additional Motivation for PACS/RIS

- Acquire, distribute and archive images from all medical imaging modalities in digital form;
- Use fast networking for transferring and accessing any part of a patient's file from any location within NIH/outside;
- Reduce the time needed for an examination to become available to the attending and referring physicians;
- Support to interhospital data exchange (teleradiology/telemedicine).
- Improve medical image interpretation with the use of image processing tools (e.g., multimodality image registration, segmentation, enhancement, etc.);
- Facilitate short- and long-term image and medical data archiving.

NIH Imaging Sciences: Imaging Modalities



PACS Viewing Stations

CR/DR/Specials/Consulting 4-hdr 2.5K AutoRads (Cemax)



MR/CT/Neuro 4-hdr 1.5K AutoRads (Cemax)



Clinical Consulting Flat-Screen 4-hdr 1.5K AutoRads (Cemax)




Neuro/ICU/Surgery AutoRads 2-hdr 1.5K (Cemax)



QC/Film Library/Clinic(s) Clinical AutoRads



Nuc Med/PET (MedImage)



MR/CT/Neuro AWWs (GE)



MR/CT/Neuro 2-hdr AWWs (GE)



Ultrasound (Acuson)



**PACS/RIS
Data Center**

**Digital Imaging
Network**

Dry Printers



Color Laser



Wet Printers (Kodak 2180s)



Web Servers



**200 comm ports
for 2,000 users**



**Clinical
Access
Server**



**50 clinical
users**




Any validated DICOM Q/R




Goldszal et al., 2000


Nuc Med/PET Mini-PACS



Ultrasound Mini-PACS




PACS/RIS Data Management Center




Digital Imaging Network (DIN)

- Redundant (at highest level) meshed network;
- High-performance, gigabit and multi-layer intelligence (up to 150Mpps);
- 4GB backbone (full-duplex), 32GB backplane;
- Support multiple protocol and VLAN switching;
- Initially, network with 225 copper drops; each drop consists of 2 CAT5 cables using redundant wire paths;
- Fiber cabling for 25 drops;
- Cisco Works 2000 network monitoring package;
- 24x7 coverage by Kodak;
- Two Pix firewall servers (full redundancy); encryption soft; VPN soft; Secure ID; PACS/RIS security policy


PACS: Web Solution



PACS: Support for Image Processing



PACS/RIS/HIS Workflow: Enterprise-wide Integration



Radiology Information System

Purpose:

To connect the appropriate persons, knowledge, and resources at the appropriate time and location to achieve the optimal health outcome.

What Capabilities are Available to CC/NIH?

- Diagnostic Radiology
- Radiology Tracking
- Departmental Scheduling
- Discern Explorer and Discern Expert
- Voice Recognition
- Interfaces
- PACS integration

Diagnostic Radiology Components

Components of Diagnostic Radiology to help manage its services.

- ◆ **Departmental Patient Registration**
- ◆ **Departmental Order Entry and Inquiry**
- ◆ **Image Management**
- ◆ **Report Processing**
- ◆ **Radiologist Desktop**
- ◆ **Quality Management**
- ◆ **Mammography QA and Patient Follow-Up**
- ◆ **Workload and Management Reporting**
- ◆ **Bar Code Reading capability**
- ◆ **Departmental Charting**
- ◆ **Charge Capture/Reporting**
- ◆ **System Operations**
- ◆ **Data Archive Management**
- ◆ **System Support**
- ◆ **Database Administration**
- ◆ **Accession Linking**

Radiology Tracking

Maintains information about the location of patients within the radiology service. Radiology Tracking allows a radiology service to keep a record of the following activities:

- The time transport was called to transfer a patient
- The time transport left to transfer a patient
- The time a patient arrived
- The time an exam was started and completed
- The time a patient entered or exited a location
- The time a patient left the radiology service

Departmental Scheduling

- » Appointment requests can be scheduled manually via an online electronic “appointment book” or by reviewing and modifying a suggested appointment schedule created by the system based on procedure and resource constraint databases.
- » Appointments can be scheduled for up to 23 months into the future.
- » For each scheduled procedure, the system can track a reason for performing the procedure, default comment templates, a notification message, special instructions, a transport code, a referring location, and free-text appointment comments.
- » User-defined, procedure-specific instructions for patient preparation can be displayed at the time an appointment is booked and can be printed on the patient’s schedule.
- » One-time, series, and dependent multiple-step appointment booking is supported.
- » Printing of appointment reminders, cancellation notices, reschedule notices, and missed appointment notices in a user-defined format is supported.
- » When a non-Cerner system is providing requests for scheduled procedures through an orders interface, Departmental Scheduling Management optionally can notify that system when those procedures or orders are scheduled, rescheduled, or canceled, via Appointment Notifications Outgoing.

Discern Explorer/Expert

Explorer

- **Retrospective decision support**
- **4GL, SQL-like language**
- **Access via command line or graphical user-friendly front end tool**

Expert

- **Prospective decision support using “if, and, or, then” rules**
- **Clinical/operational/administrative use**

Discern Expert Example

If

Carotid Angiogram is Ordered

And

Allergy to Iodine is Present

Or

Creatinine and BUN are High

Or

Taking Heparin or Coumadin

Or

Pro-Times not Ordered or Elongated

Then


Notify Physician of Possible Adverse Effects Due to Decreased Renal Function or Allergy to Contrast Media



Radiology Management Reports

- Turn-around time for completion of exams
- Room and equipment utilization rates
- Exam frequencies by NIH protocol number
- Workload by employee
- Patient waiting time
- Turn-around time for diagnostic reports
- **Over 500 other detailed and customizable reports are available from the RIS**

Typical RIS Architecture



Benefits/Opportunities Provided by the RIS

(Improvements in Quality, Data Access and/or Cost Reductions)

- ◆ Better patient follow up via Mammography QA and Patient Follow-Up
- ◆ Better resource (people, rooms, machines) scheduling via conflict checking
- ◆ Reduced waste and repeat films
- ◆ Reduced incidence of misplaced and overdue films (Film Management)
- ◆ Reduced booking and sequencing errors through use of "Overbook" and "Conflict" warnings
- ◆ Reduced turnaround times (reduced patient waits)
- ◆ Near 100% capture of orders and orderables
- ◆ Improved workload and management reporting that provides the factual basis for corrective actions
- ◆ Patient care improvements obtained by use of duplicate checking, Discern Expert rules and clinical alerts
- ◆ Seamless integration with PACS